**Prótese Valvular Biológica  
INOVARE® SAFESYNC –  
B-Delivery System**

**INOVARE® SAFESYNC  
Biological Valvular Prosthesis –  
B-Delivery System**

**Prótesis Valvular Biológica  
INOVARE® SAFESYNC –  
B-Delivery System**

Cod. Manual: 623076 - D-1

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## 1. PRODUCT DESCRIPTION

The INOVARE® SAFESYNC Biological Valvular Prosthesis – B-Delivery System comprises the following components:

- INOVARE® SAFESYNC Biological Valvular Prosthesis
- B-Delivery System
- Pre-Dilation Balloon Catheter
- Aureus® Transfemoral Expandable Introducer
- Crimper

The medical devices listed above can be used as replacement components for exclusive use in the system.

### 1.1. INOVARE® SAFESYNC Biological Valvular Prosthesis

The INOVARE® SAFESYNC Biological Valvular Prosthesis (Figure 1) is a permanent implant device, classified as a balloon-expandable transcatheter biological cardiac valve prosthesis. It is intended for transcatheter aortic valve replacement (TAVR).

Figure 1: INOVARE® SAFESYNC Biological Valvular Prosthesis as observed (from left to right) in isometric, frontal and superior view.

For the implantation, the valve prosthesis is compressed over the inflation balloon of its delivery system, and, loaded on it, is navigated through the anatomy of the patient and delivered to the intended implant site, where the balloon is inflated to expand and release the valve prosthesis.

The component is presented in six different models, whose dimensions are presented in Table 1.

Table 1: INOVARE® SAFESYNC Biological Valvular Prosthesis models and their respective specifications.

Model	External Diameter (mm)	Internal Diameter (mm)
20	20	18
22	22	20
24	24	22
26	26	24
28	28	26
30	30	28

The sizing recommendations presented in Table 2 are based on the diameter of the native aortic valve annulus, as measured by computed tomography (CT). The anatomical factors of the patient and multiple imaging modalities should be considered during the selection of the diameter of the valve prosthesis. Moreover, risks associated with undersizing and oversizing should be evaluated.

Table 2: Sizing recommendations for the valve prosthesis.

Native aortic valve annulus (CT)		Dimensions of the recommended valve prosthesis	
Area (mm <sup>2</sup> )	Area derived diameter (mm)	Outer diameter (mm)	Area (mm <sup>2</sup> )
251 – 299	17,9 – 19,5	20	314
299 – 362	19,5 – 21,5	22	380
362 – 431	21,5 – 23,4	24	452
425 – 506	23,3 – 25,4	26	531
493 – 586	25,0 – 27,3	28	616
565 – 673	26,8 – 29,3	30	707

The valve prosthesis is composed of an expandable metallic structure made of cobalt-chromium alloy, onto which the polyester fabric skirt and bovine pericardium cusps are sewn with polyester sutures. Tantalum radiopaque markers are embedded in the expandable metal structure.

### 1.1.1. Expandable metal structure

The expandable metal structure is made of cobalt-chromium alloy conformed into a ring-shaped geometry, consisting of successive V's (Figure 2). It is produced by laser cutting and finished by electropolishing. The polyester skirt and the bovine pericardium cusps are sewn onto this structure using polyester suture. Tantalum radiopaque markers are embedded in the three poles that are distributed radially equidistant along the metal structure.


Figure 2: Ring-shaped expandable metal structure with poles distributed radially equidistant along it.

### 1.1.2. Radiopaque markers

Radiopaque markers are microspheres of the chemical element Tantalum, highlighted in yellow in Figure 3. The radiopacity of the markers provides fluoroscopic guidance during the implant, allowing the necessary traceability for the safe and accurate release of the INOVARE® SAFESYNC Biological Valvular Prosthesis.

Figure 3: Radiopaque markers embedded in a post of the metal structure (highlighted in black).

### 1.1.3. Polyester fabric skirt and polyester sutures

The skirt (Figure 4) is made of polyester, a material with low permeability, allowing greater sealing in the aortic ring and reduction of paravalvular leaks (PVL). The sutures used to sew the skirt and the cusps onto the metal structure are also made of polyester, due to the excellent mechanical resistance of the material, allowing a robust union between metallic structure, skirt and cusps.

Figure 4: Polyester skirt (white).

### 1.1.4. Bovine pericardium cusps

The bovine pericardium used in making the cusps is obtained immediately after the animal is slaughtered in slaughterhouses previously examined by the Federal Inspection Service (SIF). Braile Biomédica's bovine pericardium is free of transmissible spongiform encephalopathies and foot-and-mouth disease.

The material is collected from animals with a maximum age of 36 months, and only pericardium rich in collagen fibers and free of fat are selected for the manufacture of the valve prosthesis.

The material collected undergoes chemical treatment and is placed in a conservative solution of 4% formaldehyde to be subsequently submitted to quality control and released for the manufacture of the valve prosthesis.

The bovine pericardium, after its treatment and release, is sutured with polyester threads to the metallic cobalt-chromium support previously covered with polyester fabric, thus forming the valve cusps (Figure 5).


Figure 5: Details of the bovine pericardium cusps a) as seen from the aortic face (top) and b) from the ventricular face (bottom).

## 1.2. B-Delivery System

The B-Delivery System is used for the delivery and implantation of the INOVARE® SAFESYNC Biological Valvular Prosthesis in TAVR procedures with transfemoral access. It consists of an inflation balloon mounted in a deflectable sheath, which is coupled to a handle that actuates the delivery system's mechanisms. The valve prosthesis is compressed over the inflation balloon and loaded into the deflectable sheath, and the delivery system is navigated through the patient's anatomy to the intended implant site, where the balloon is inflated to deploy the valve prosthesis.

Figure 6 presents the delivery system and its parts in detail.


Figure 6: B-Delivery System detailing its parts.

Chart 1 details the functions of each of the delivery system's parts.

Chart 1: B-Delivery System parts and their respective function.

Delivery System Parts	Function
1. Inflation balloon	Deployment of the valve prosthesis
1.1. Atraumatic tip	Prevention of vascular damage during delivery
1.2. Radiopaque markers	Traceability and guidance for the positioning of the valve delivery system
2. Deflectable sheath	Navigation/delivery of the valve prosthesis Flexion mechanism for facilitated navigation in the aortic arch
3. Loader	Protection of the valve prosthesis during passage through the haemostatic valves of the expandable introducer, avoiding damage and loss of positioning
4. Handle	Activation of delivery system mechanisms
4.1. Flex wheel	Triggering sheath deflection
4.2. Flush port	Flushing of the deflectable sheath's lumen
4.3. Haemostatic valves	Prevention of blood leakage
4.4. Positioning rod unlocking button	Allow movement of the positioning rod Prevention of undesirable displacement of the inflation balloon (when locked, i.e., unpressed)
4.5. Fine adjustment wheel	Allow the inflation balloon to be retrieved into the deflectable sheath, providing refinement of the positioning of the valve prosthesis for deployment
5. Positioning rod	Allow adjustment of the inflation balloon's position in relation to the deflectable sheath
5.1 Hub	Guidewire lumen port (straight port) Inflation balloon lumen port (side port)

The delivery system is presented in six different models, each corresponding to one of the six valve prosthesis models available. Its dimensions, characteristics of balloon inflation and compatibility with valve prosthesis models are presented in Table 3. All delivery system's models are designed to be compatible with 0.035" (0.889 mm) guidewires.

Table 3: Dimensional and inflation specifications of the delivery system and compatibility with the valve prosthesis.

Delivery system model	Balloon diameter x Useful balloon length x Total device length	Nominal Volume (NV)	Nominal Pressure (NP)	Rated Burst Pressure (RBP)	Compatible valve prosthesis model
20	20 mm x 3.5 cm x 110 cm	11 mL	4 atm	8 atm	20
22	22 mm x 3.5 cm x 110 cm	15 mL	4 atm	8 atm	22
24	24 mm x 3.5 cm x 110 cm	21 mL	4 atm	8 atm	24
26	26 mm x 3.5 cm x 110 cm	24 mL	4 atm	8 atm	26
28	28 mm x 3.5 cm x 110 cm	27 mL	4 atm	8 atm	28
30	30 mm x 3.5 cm x 110 cm	30 mL	4 atm	5 atm	30

The recommendations for vascular access of the delivery system with the loaded valve prosthesis are presented in Table 4. Due to the dimensional differences between the models of the valve prosthesis, its dimensions when loaded in the delivery system also diverge, requiring specific arterial calibers for safe vascular access. In addition to the dimensions, anatomical factors of the patient such as tortuosity, presence of obstructions and calcification should be considered as they impact access safety directly.

Table 4: Compatibility between product model and vascular access diameter.

Model	Vascular Access Diameter (mm)
20	≥ 6,50
22	≥ 6,50
24	≥ 6,50
26	≥ 6,50
28	≥ 7,00
30	≥ 7,00

### 1.3. Pre-Dilation Balloon Catheter

The Pre-Dilation Balloon Catheter is intended for expanding the patient's stenosed aortic valve before valve prosthesis implantation (i.e., pre-dilation), which results in calcium fracture and separation of fused cusps, increasing the orifice area to allow for easier access of the delivery system with the loaded valve. It also allows for greater uniformity in the expansion of the valve prosthesis during the implant by reducing radial counterforces, which, in turn, refines the positioning of the implant and reduces PVL.

Figure 7 shows the Pre-Dilation Balloon Catheter and its parts.


Figure 7: Pre-Dilation Balloon Catheter and its parts.

The pre-dilation balloon is presented in diameters of 18 and 20 mm, being compatible with all models of the INOVARE® SAFESYNC Biological Valvular Prosthesis (the choice of model is based on the patient's anatomical conditions). Table 5 shows the dimensional and inflation specifications of the pre-dilation balloon.

Table 5: Dimensional and inflation specifications of the pre-dilation balloon.

Model	Balloon diameter x Useful balloon length x Total device length	Nominal Volume (NV)	Nominal Pressure (NP)	Rated Burst Pressure (RBP)
18	18 mm x 4.0 cm x 110 cm	11 mL	6 atm	8 atm
20	20 mm x 4.0 cm x 110 cm	14 mL	5 atm	8 atm

#### 1.4. Aureus® Transfemoral Expandable Introducer

The Aureus® Transfemoral Expandable Introducer (Anvisa Register No. 10159030110) consists of a hydrophilic expandable sheath, a dilator to facilitate transfemoral percutaneous vascular access and an expander to promote radial expansion of the sheath after access is established, as illustrated in Figure 8.

Figure 8: Aureus® Transfemoral Expandable Introducer: (1) expandable sheath, (2) dilator and (3) expander.

Figure 9 shows the parts of the Aureus® Transfemoral Expandable Introducer. It illustrates the device in its percutaneous vascular access configuration, with the dilator fully inserted in the lumen of the expandable sheath. This configuration presents a vascular access diameter of 14 Fr. After access is established, the dilator is removed and the expander is inserted, promoting the radial expansion of the sheath. The expander can then be removed, and now the expandable sheath allows the passage of devices of up to 22 Fr.


Figure 9: Description of the components of the Aureus® Transfemoral Expandable Introducer.

The dimensional specifications of the Aureus® Transfemoral Expandable Introducer are presented in Table 6.

Table 6: Dimensional specifications of the Aureus® Transfemoral Expandable Introducer.

Internal Diameter		External Diameter		Useful length
Expandable Sheath				
Fr	mm	Fr	mm	mm
14	4,67	19	6,33	260
Dilator				
Pol. (")	mm	Fr	mm	mm
0,035	0,89	14	4,67	493
Expander				
Pol. (")	mm	Fr	mm	mm
0,035	0,89	22,5	7,50	353

## 1.5. Crimper

The Crimper is used to radially compress (i.e., crimp) the valve prosthesis on the inflation balloon of the delivery system. The crimping process allows the valve prosthesis to have its diameter reduced and be firmly attached to the inflation balloon, so that it can be safely delivered to the intended implant site. Its parts are detailed in Figure 10.


Figure 10: Crimper and its parts.

The Crimper's operating principle is based on the compression caused by the movement of the handle, which leads to the closure of its compression area, thus triggering the crimping mechanism (Figure 11).


Figure 11: Crimper operation: a) rotation of the handle leading to the closure of the compression area, and b) compression of the valve prosthesis on the inflation balloon of the delivery system.

## 2. HOW SUPPLIED

The INOVARE® SAFESYNC Biological Valvular Prosthesis – B-Delivery System is available in the models presented in Table 7, which details the main dimensions of each component of the system as well as its intercompatibility.

Table 7: Components of the INOVARE® SAFESYNC Biological Valvular Prosthesis - B-Delivery System and compatibility between their model.

Component	System Model 20	System Model 22	System Model 24	System Model 26	System Model 28	System Model 30
INOVARE® SAFESYNC Biological Valvular Prosthesis	20 mm	22 mm	24 mm	26 mm	28 mm	30 mm
B-Delivery System	20 mm	22 mm	24 mm	26 mm	28 mm	30 mm
Pre-Dilation Balloon Catheter	18 or 20 mm					
Aureus® Transfemoral Expandable Introducer	14 Fr x 260 cm					
Crimper	-					

The components listed in Table 7, inside their respective individual packaging, are packed in two corrugated carton boxes, in the following arrange:

- one corrugated carton box contains (depending on the System's model):
  - ✓ 01 (one) INOVARE® SAFESYNC Biological Valvular Prosthesis model 20, 22, 24, 26, 28 OR 30 mm
  - ✓ 01 (one) Pre-Dilation Balloon Catheter model 18 OR 20 mm
  - ✓ 01 (one) Aureus® Transfemoral Expandable Introducer
  - ✓ 01 (one) Crimper
- the other corrugated carton box contains (once again depending on the System's model):
  - ✓ 01 (one) B-Delivery System model 20, 22, 24, 26, 28 OR 30 mm (corresponding to the model of the valve prosthesis that it accompanies).

Sections 2.1through 2.5 detail the individual packaging of each component.

## 2.1. INOVARE® SAFESYNC Biological Valvular Prosthesis

The valve prosthesis is supplied sterile (sterilized by chemical sterilization) and non-pyrogenic, mounted on a support structure and packed in a polycarbonate container with 4% formaldehyde solution. The container is accommodated in an expanded polystyrene (EPS - Styrofoam) cradle, and packed in a cardboard paper box. Two temperature indicators are coupled to the Styrofoam cradle, and the cardboard paper box has a transparent plastic covered side opening that allows the visualization of the indicators. In the cardboard paper box, one Instructions Manual, one implant file and six traceability labels are packed with the valve prosthesis.

The packaging is presented in detail in Figure 12.


Figure 12: Packaging of the INOVARE® SAFESYNC Biological Valvular Prosthesis: a) valve prosthesis mounted on its support assembly inside polycarbonate container, b) Styrofoam cradle with temperature indicators, and c) cardboard box with side opening covered in transparent plastic for visualization of temperature indicators.

## 2.2. B-Delivery System

The B-Delivery System is supplied sterile (sterilized by ethylene oxide - ETO) and non-pyrogenic. The component is accommodated in a polyethylene cradle, with protectors for the guidewire lumen and for the deflectable sheath, packed in a double microbial barrier polyethylene envelope, and finally packed in a cardboard paper box.

The packaging is presented in detail in Figure 13.


Figure 13: Packaging of the B-Delivery System: a) delivery system, b) delivery system with protectors for the lumen of the guide wire and for the deflectable sheath, c) accommodated in a polyethylene cradle, d) double microbial barrier envelope, and e) cardboard paper box.

### **2.3. Pre-Dilation Balloon Catheter**

The Pre-Dilation Balloon Catheter is supplied sterile (sterilized by ETO) and non-pyrogenic. It is arranged in a polyethylene tray with vinyl acetate (EVA) holders, packed in a microbial barrier polyethylene envelope, and then packed in a cardboard paper box. Its packaging is presented in detail in Figure 14.


Figure 14: Packaging of the Pre-Dilation Balloon Catheter: a) Pre-Dilation Balloon Catheter, b) arranged in a tray with EVA supports, c) microbial barrier envelope, d) cardboard paper box.

### **2.4. Aureus® Transfemoral Expandable Introducer**

The Aureus® Transfemoral Expandable Introducer is supplied sterile (sterilized by ETO) and non-pyrogenic. It is arranged in a polyethylene tray with EVA holders, packed in a microbial barrier polyethylene envelope, and then packed in a cardboard paper box. Its packaging is presented in detail in Figure 15.


Figure 15: Packaging of the Aureus® Transfemoral Expandable Introducer: a) Aureus® Transfemoral Expandable Introducer, b) arranged in a tray with EVA supports, c) microbial barrier envelope, d) cardboard paper box.

## 2.5. Crimper

The Crimper is supplied sterile (sterilized by ETO) and non-pyrogenic. It is packed in a double microbial barrier polyethylene envelope and then packed in a corrugated carton box.

Its packaging is shown in detail in Figure 16.


Figure 16:  
Crimper packaging: a) Crimper,  
b) double microbial barrier envelope,  
c) corrugated carton box.

### 3. INTENDED USE

The INOVARE® SAFESYNC Biological Valvular Prosthesis – B-Delivery System is indicated for the treatment of symptomatic patients with severe aortic valve stenosis (AoS) who, evaluated by a multidisciplinary team of physicians specializing in valve diseases, have a predicted risk of mortality from conventional surgery  $\geq 8\%$  at 30 days, based on the Society of Thoracic Surgeons risk score (STS score) and other comorbidities not scored by the score, being judged appropriate for the treatment of transcatheter aortic valve replacement (TAVR) with access transfemoral.

### 4. CONTRAINDICATIONS

The INOVARE® SAFESYNC Biological Valvular Prosthesis – B-Delivery System is contraindicated in patients who cannot tolerate an anticoagulation/antiplatelet regimen or who have active bacterial endocarditis or other active infections.

## 5. WARNINGS AND PRECAUTIONS

### 5.1. Warnings

- Read the instructions manual before use.
- Product intended for single use only. Do not reuse and do not reprocess.
- The medical devices that compose the INOVARE® SAFESYNC Biological Valvular Prosthesis – B-Delivery System are supplied sterile. Do not resterilize the devices by any method. There is no data to support the sterility, non-pyrogenicity and functionality of the devices after resterilization.
- Prior to the procedure, the INOVARE® SAFESYNC Biological Valvular Prosthesis and the B-Delivery System should be examined to ensure that their sizes are compatible and suitable for their intended use.
- The compatibility between INOVARE® SAFESYNC Biological Valvular Prosthesis models and B-Delivery System models is essential to minimize risks and complications to the patient.
- Do not use the B-Delivery System for delivery of a valve prosthesis other than the INOVARE® SAFESYNC Biological Valvular Prosthesis. There is no data to support the compatibility of the delivery system with other valve prostheses.
- Accelerated deterioration of the valve prosthesis may occur in patients with an altered calcium metabolism.
- The correct choice of the valve prosthesis model is essential to minimize the risk of paravalvular leak, embolization, migration, residual gradient (prosthesis-patient mismatch) and/or annular rupture.
- The physician should check the correct orientation of the valve prosthesis before implantation: the aortic face should be pointing to the handle of the delivery system.
- Do not expose the valve prosthesis to extreme temperatures. Prolonged exposure to heat (30°C and higher) or cold (10°C and lower) can irreversibly damage the biological tissue. Valve prostheses exposed to extreme temperatures, or under suspicion of exposure should not be used. Observe the temperature indicators located on the side of the package to check the exposure of the valve prosthesis to temperatures beyond its specified limits. If there is exposure, the indicators will turn red. In this case, replace the valve prosthesis.

- Carefully inspect all packaging. Do not use if they are damaged, are outside their expiration date, or if the sterile barrier is compromised. Discard such devices and replace with identical ones. Specifically for the valve prosthesis, damage include cracks in the container or lid, leakage of conservation solution, and violation or lack of plastic casing and/or safety seal.
- Do not use the valve prosthesis if it is not completely submerged in the conservation solution. In this case, replace it and inform a representative of Braile Biomédica.
- The outer surface of the container should not be considered sterile. Never put it on sterile fields.
- Do not expose the valve prosthesis to solutions other than the conservation solution (formaldehyde) in which it is supplied, or the sterile physiological saline used during the rinsing process.
- Do not add antibiotics to the conservation and rinsing solutions. Do not utilize antibiotics in the valve prosthesis. The valve prosthesis is sterilized by glutaraldehyde and preserved in formaldehyde, requiring no antibiotic or fungicide.
- Do not allow the valve prosthesis to dry. It should be hydrated at all times between removal of the storage solution and delivery.
- Do not use unprotected tweezers or sharp instruments while handling the valve prosthesis. Avoid touching the biological tissue.
- Formaldehyde may cause irritation of the skin, eyes, nose and throat. Avoid prolonged or repeated exposure to, or breathing of, the solution. Use only with adequate ventilation. If skin contact occurs, immediately flush the affected area with water; in the event of contact with eyes, seek immediate medical attention.
- Do not exceed the balloon working pressure of the delivery system during inflation, as this action may result in balloon rupture, increasing the overall risk of the valve prosthesis implantation procedure.
- When the delivery system is in contact with the vascular system, it should be manipulated under high-quality fluoroscopic observation. Do not advance or retreat the delivery system unless the balloon is fully deflated by vacuum action (zero pressure). When encountering resistance during manipulation, determine the cause before proceeding.
- Do not use in the case of tampered or damaged devices.
- Do not use if the expiration date has elapsed.

- Do not use if the labelling is incomplete or illegible.
- Do not handle the devices improperly, as they may be damaged prior to use.
- After use, dispose of medical devices and packaging in accordance with hospital, administrative and/or local government policy.
- The use of excess contrast solution can lead to acute kidney injury. Monitor the patient's creatinine level before the procedure. The administration of contrast should be monitored throughout the procedure.
- Monitor the pacemaker electrode throughout the procedure to avoid the potential risk of right ventricle perforation.
- There may be a higher risk of stroke in TAVR procedures compared to balloon aortic valvuloplasty or other standard procedures for high or extreme risk patients.
- Precautions should be taken when implanting the valve prosthesis in patients with severe coronary artery disease.
- Lesions to the patient's vascular system may occur if the delivery system is deflected during removal. Completely undo deflection before removal.
- The procedure should be conducted under fluoroscopic guidance. Some procedures that employ fluoroscopy are associated with the risk of skin injury caused by radiation. These lesions may be painful, disfiguring, and long-lasting.
- Patients with the implanted valve prosthesis should be kept under anticoagulation/antiplatelet therapy, except when contraindicated by their doctor. This device has not been tested for use without anticoagulation.

## 5.2. Precautions

- Read all warnings, precautions and instructions for use carefully before using the product.
- It is necessary to make sure of the integrity of the product before its use.
- Do not use the product if it or the packaging is damaged.
- Do not use the valve prosthesis if the seal is tampered with or if the conservative solution leaks.
- Do not use the product in case of doubt.
- Always use aseptic techniques when handling of the product.

- Do not let the valve prosthesis tissue dry. Use only sterile isotonic physiological saline, immediately after removal of the valve prosthesis from the container until the moment of implantation.
- The long-term durability of the valve prosthesis has not been established. Routine consultations are advised to evaluate its performance.
- Do not aspirate the steam from the formaldehyde solution used for conservation.
- Take care when using cutting instruments as it may damage the biological tissue of the valve prosthesis.
- The outer side of the valve prosthesis' container should not be considered sterile. Never put it on the instrumentation table.
- Always maintain a backup unit of each component during the procedure.
- The physician is responsible for the final decision on the product's model.
- The prosthesis implantation should be done according to the surgeon's technique.
- Store the product in a clean, dry, weather-free environment at temperatures above 10°C and below 30°C. Never freeze or expose to temperatures above 30°C. Freezing or prolonged exposure to heat can irreversibly damage biological tissue. Products that are subject to extreme temperatures, or under suspicion of exposure, should not be used.
- Check the expiration date before using the product. The expiration date indicated on the product label must be strictly respected.
- Product intended for single use only. It should not be resterilized or reused, as it may compromise its performance, in addition to increasing the risk of improper sterilization with possible contamination.
- The sale and use of this product are restricted to physicians only.
- During the surgical procedure it is recommended that the patient be kept on anticoagulation/antiplatelet therapy, according to the medical prescription, to minimize the risk of thrombosis or thromboembolic occurrences.
- Do not overinflate the delivery system's balloon during the implantation of the valve prosthesis, as this may prevent proper coaptation of leaflets, affecting the valve prosthesis' performance.
- The delivery system should be used only by physicians with experience in vascular procedures for implantation of transcatheter aortic valve prostheses with femoral percutaneous access.

The safety and effectiveness of the INOVARE® SAFESYNC Biological Valvular Prosthesis - B-Delivery System were not evaluated for patients with the following characteristics/comorbidities:

- Recent transient stroke or ischemic attack (over the previous 6 months)
- Anemia (hemoglobin < 10 g/dl), thrombocytopenia (< 100,000 cel/mm<sup>3</sup>) or hyperthrombocytosis (> 700,000 cel/mm<sup>3</sup>);
- Absence of calcification in the native valve;
- Valve ring area-derived diameter smaller than 17.3 or greater than 28.6 mm;
- Native valvar dysfunction of mixed component (stenosis and insufficiency) with predominant regurgitation or ≥ 3+;
- Coronary disease requiring elective revascularization during or after the valve implantation procedure;
- Active peptic disease, gastrointestinal bleeding < 3 months or previously diagnosed hemorrhagic diasis;
- Endocarditis < 12 months;
- Evidence of myocardial infarction in an interval shorter than one month;
- Evidence of intracardiac mass (tumor, thrombus or vegetation);
- Life expectancy less than 12 months due to non-cardiac disease or other comorbidities;
- Left ventricle ejection fraction < 30%;
- Pregnancy;
- Hypersensitivity to cobalt, chromium, nickel, molybdenum, titanium, manganese and polymeric materials;
- Hypersensitivity or contraindication to the administration of heparin, ticlopidine or clopidogrel;
- Hypersensitivity or contraindication to radiopaque contrast medium;
- Hemodynamic instability requiring vasoactive drugs or circulatory support;
- Marked aortic, mitral or tricuspid insufficiency;
- Acute renal failure with serum creatinine that has not yet returned to baseline levels;
- Chronic dialysis renal failure or with serum creatinine levels > 3.0 mg/dL (265 µmol/L);
- Obstructive hypertrophic cardiomyopathy or severe obstruction of the left ventricular outflow tract, with no possibility of balloon dilation;

- Need for chronic anticoagulation for other causes;
- Severe obstruction of the outflow route, without possibility of balloon dilation;
- Presence of sessile and institable atheromas in the ascending aorta and/or aortic arch detected by imaging methods;
- Previous aortic or mitral valve procedure (surgical or catheter, excluding balloon aortic valvuloplasty);
- Valve procedure with clinical need to perform in urgent or emergency time (non-elective);
- Clinical or biological signs of infection with systemic repercussion;
- Noncalcified native valve;
- Iliac-femoral access route with extreme tortuosity or calcification or reduced luminal diameter (typically < 7 mm), which prevents the safe progression of the arterial introducer and delivery system, according to the judgment of the responsible physician.

### 5.3. MR SAFETY AND COMPATIBILITY

MR Conditional

Non-clinical testing has shown that the INOVARE® SAFESYNC Biological Valvular Prosthesis is conditional for Magnetic Resonance Imaging (MRI), which means that its safety has been demonstrated for a MRI environment under the following conditions:

- static magnetic field of 3.0 T;
- spatial gradient field of 22900 gauss/cm (229 T/m, extrapolated value) or less;
- product of static magnetic field and spatial gradient (force product) of 435 T<sup>2</sup>/m (extrapolated value) or less;
- maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2.4 W/kg for 15 minutes of scanning in Normal Operation Mode.
- Under the scan conditions defined above, the INOVARE® SAFESYNC Biological Valvular Prosthesis is expected to produce a maximum temperature rise of 0.3 °C after 15 minutes of continuous scanning.
- In non-clinical testing, the image artifact caused by a worst-case implant extends as far as 21.2 mm from the implant for gradient echo images when scanned in a 3.0 Tesla MRI system.

The valve prosthesis has not been evaluated in MRI systems other than 3.0 T.

## 6. POTENTIAL ADVERSE EVENTS

Complications may occur peri- or post-procedure. Procedures requiring treatment should be performed only by physicians with due knowledge regarding the possible associated complications, which include, but are not limited to:

### 6.1. Events Related to Clinical Aspects or the Procedure

- Cardiovascular complications including:
  - ✓ Amputation;
  - ✓ Aneurysm;
  - ✓ Angina;
  - ✓ Arrhythmias, including ventricular fibrillation (VF) and ventricular tachycardia (VT);
  - ✓ Arterial and venous occlusion;
  - ✓ Arteriovenous fistula (AV) or pseudoaneurysm;
  - ✓ Cardiogenic shock;
  - ✓ Conversion to conventional surgery;
  - ✓ Dissection, perforation or vascular ruptures;
  - ✓ Heart failure / infarction;
  - ✓ Hematoma;
  - ✓ Hypertension or hypotension;
  - ✓ Occlusion of collateral vessels;
  - ✓ Perforation or rupture of cardiac structures;
  - ✓ Pericardial effusion or cardiac tamponade;
  - ✓ Peripheral ischemia or nerve damage;
  - ✓ Syncope;
  - ✓ Thromboembolism;
  - ✓ Thrombosis;
  - ✓ Tissue necrosis;
  - ✓ Vessel trauma needing repair or surgical intervention;
  - ✓ Vascular ischemia;

- ✓ Vasospasm;
- ✓ Vasovagal response.
- :Neurological complications, including:
  - ✓ Blindness;
  - ✓ Neuropathy;
  - ✓ Paralysis;
  - ✓ Paraparesis;
  - ✓ Paresthesia;
  - ✓ Stroke;
  - ✓ Transient ischemic accidents.
- Pulmonary complications, including:
  - ✓ Embolism, distal (gas embolisms, tissues or thrombotics);
  - ✓ Pulmonary edema;
  - ✓ Respiratory failure or respiratory arrest.
- Allergic reaction to antithrombotic therapy, contrast medium or anesthesia;
- Anemia;
- Bleeding;
- Compartment syndrome;
- Complications related to wound healing;
- Death;
- Edema;
- Infection and fever;
- Inflammation;
- Pain or changes in the access site;
- Renal failure.

## 6.2. Events Related to the Device

- Allergic/immunological reaction to the implant;
- Atrial fibrillation, atrial flutter;
- Bleeding requiring transfusion or intervention;
- Cardiac arrest;

- Cardiogenic shock;
- Compression of the conduction system including left atrioventricular block, which may require a provisional pacemaker;
- Coronary occlusion;
- Difficulties in the introduction, navigation and removal of the delivery system;
- Dissection, rupture, trauma of the aortic ring and surrounding structures including the ascending aorta, coronary ostium, and ventricular septum;
- Emergency cardiac surgery;
- Explant of the valve prosthesis;
- Fracture of the stent of the valve prosthesis;
- Heart failure or reduced cardiac output;
- Hemolysis;
- Implantation of the valve prosthesis in an unintended location;
- Inability to implant the valve prosthesis;
- Inaccurate implant (e.g., inadequate guidance of the valve prosthesis);
- Infection, fever, septicemia, abscess, endocarditis;
- Mechanical failure of the delivery system, including balloon rupture and tip separation;
- Mediastinal bleeding;
- Mediastinitis;
- Migration of the valve prosthesis, inadequate positioning or embolization that require intervention;
- Mitral valve injury;
- Paravalvar leak;
- Regurgitation of the valve, paravalvular or transvalvular;
- Rupture of graft material;
- Silent cerebral ischemia, stroke, transient ischemic accident, cognitive impairments;
- Structural deterioration of the valve prosthesis (wear, fracture, calcification);
- Thrombosis of the valve prosthesis.

## **7. PACKAGING, STORAGE AND TRANSPORTATION**

The components of the INOVARE® SAFESYNC Biological Valvular Prosthesis - B-Delivery System must be preserved in their respective individual packaging, and packed in the respective cardboard boxes, as supplied. The product must be stored in a clean, dry, weather-free environment, at temperatures above 10 °C and below 30 °C.

Storage also requires proper care to prevent falls and possible damage to the product. For transport, it should never be cooled, frozen or exposed to temperatures above 30°C and below 10°C, and exposure to water, direct sunlight, extreme temperatures and high humidity should be avoided.

## **8. INSTRUCTIONS FOR USE**

### **A. Minimum Requirements For Product Use**

- A.1. Preoperative exams to assess patient eligibility include:
  - A.1.1. Complete blood count
  - A.1.2. Serum sodium dosage
  - A.1.3. Serum potassium dosage
  - A.1.4. Serum ionic calcium dosage
  - A.1.5. Serum urea dosage
  - A.1.6. Serum creatinine dosage
  - A.1.7. Fasting blood glucose
  - A.1.8. Serum magnesium dosage
  - A.1.9. Serum phosphorus dosage
  - A.1.10. Prothrombin time
  - A.1.11. Activated partial thromboplastin time (APTT)
  - A.1.12. Transthoracic echocardiogram (TTE)
  - A.1.13. Transesophageal echocardiogram (TEE)
  - A.1.14. Coronary angiography
  - A.1.15. Doppler ultrasound of the iliofemoral system
  - A.1.16. Doppler ultrasound of the carotid and vertebral arteries
  - A.1.17. Chest radiography

A.2. The following information should be available for surgical planning:

A.2.1. Information obtained through echocardiogram:

- A.2.1.1. Aortic valve diameter
- A.2.1.2. Aortic valve gradient
- A.2.1.3. Degree of aortic valve insufficiency
- A.2.1.4. Mitral valve dynamics
- A.2.1.5. Ejection fraction
- A.2.1.6. Presence of thrombi
- A.2.1.7. Presence of unstable atherosclerotic plaques in the ascending aorta, aortic arch and descending aorta

A.2.2. Information obtained through Electrocardiography-triggered Computed Tomography (ECG-triggered CT):

- A.2.2.1. Native aortic valve structural visualization
- A.2.2.2. Aortic valve plan angulation;
- A.2.2.3. Native aortic valve diameter;
- A.2.2.4. Left ventricular outflow tract (LVOT) diameter;
- A.2.2.5. Membranous septum visualization;
- A.2.2.6. Height of the coronary arteries relative to the valve plane;
- A.2.2.7. Relationship between leaflets lengths and coronary arteries heights;
- A.2.2.8. Valsalva sinus diameter;
- A.2.2.9. Sinotubular junction diameter;
- A.2.2.10. Ascending aorta diameter;
- A.2.2.11. Aortic-mitral angle (should be greater than or equal to 120 degrees);
- A.2.2.12. Presence of thrombi and/or calcium plaques in the ascending aorta, aortic arch and descending aorta;
- A.2.2.13. Assessment of the vascular system in relation to iliac and femoral arteries diameter, calcification and tortuosity.

**WARNING: THERE SHOULD BE BASIC MONITORING WITH CARDIOSCOPE, OXIMETER, CAPNOGRAPHY, INVASIVE BLOOD PRESSURE, NASOPHARYNGEAL TEMPERATURE, DIURESIS AND CENTRAL VENOUS PRESSURE DURING THE SURGICAL PROCEDURE. ADHESIVE EXTERNAL DEFIBRILLATION PLATES MUST BE AVAILABLE AND POSITIONED TO ALLOW FOR A QUICK CONVERSION TO SIDE THORACOTOMY OR EMERGENCY STERNOTOMY.**

**WARNING: TRANSESOPHAGEAL AND FLUOROSCOPIC ECHOCARDIOGRAPHIC IMAGING EQUIPMENT WITH DIGITAL RESOURCES MUST BE AVAILABLE.**

**ATTENTION: MEDICAL STAFF MUST BE PREPARED FOR IMMEDIATE CONVERSION FROM TRANSCATHETER TREATMENT TO CONVENTIONAL SURGERY.**

A.2.3. Required experience in order to use the device

In order to use the INOVARE® SAFESYNC Biological Valvular Prosthesis – B-Delivery System, the implant physician must be experienced in:

- A.2.3.1. Transfemoral access and catheterization
- A.2.3.2. Balloon aortic valvuloplasty
- A.2.3.3. Implantation of aortic valve prosthesis via transfemoral approach
- A.2.3.4. Knowledge and training on the implantation of INOVARE® SAFESYNC Biological Valvular Prosthesis via transfemoral approach using the B-Delivery System

A.2.4. Required materials and equipment

- A.2.4.1. INOVARE® SAFESYNC Biological Valvular Prosthesis - B-Delivery System (Table 7)
- A.2.4.2. Cardiac catheterization lab / hybrid operating room
- A.2.4.3. Standard cardiac catherization lab materials and equipment
- A.2.4.4. Standard operating room materials and equipment for the implantation of cardiac valve prosthesis
- A.2.4.5. Fluoroscopy system appropriate for use in percutaneous coronary intervention
- A.2.4.6. Transesophageal or transthoracic echocardiography capabilities

- A.2.4.5. 18G Seldinger needle
- A.2.4.6. 0.035" x 260 cm (or longer) extra-stiff guidewire, such as Safari extra small (Ref. 39407XS H749 0), Safari small (Ref. H749 39407S 0), Confida (Ref. GWBC30), Lunderquist DC (Ref. TSCMG-35-260-E-LESDC), Amplatz (Ref. THSCF-35-260-3-AES)
- A.2.4.7. Valved introducers and diagnostic catheters commonly used in hemodynamic
- A.2.4.8. Pacemaker with rapid pace function of up to 220 bpm
- A.2.4.9. 3 liters or more of physiological saline (at room temperature)
- A.2.4.10. Five sterile rinsing basins with a minimum capacity of 500 mL/cc for the physiological saline
- A.2.4.11. 50 mL or more of non-ionic iodized radiopaque contrast solution.
- A.2.4.12. One sterile basin with a minimum capacity of 150 mL / cc for the mixture of physiological saline and radiopaque contrast solution.
- A.2.4.13. Sterile table with waterproof cover for devices preparation
- A.2.4.14. Two balloon inflation devices which graduations compatible with the nominal inflation volumes of the delivery system balloon (Table 3) and the pre-dilation balloon (Table 5), such as the VDK-12 60 mL Balloon Inflation Device (Ref. 612154)
- A.2.4.15. Two high-pressure 3-way stopcocks
- A.2.4.16. 20 mL / cc syringe
- A.2.4.17. 60 mL / cc luer-lock syringe
- A.2.4.18. Forceps
- A.2.4.19. Scalpel

Before beginning the preparation of any of the INOVARE® SAFESYNC Biological Valvular Prosthesis – B-Delivery System components, carefully inspect all packagings. Do not use if damage or defects are observed, if the expiration date has elapsed or if the sterile barrier is compromised. In that case, dispose of said devices and replace with identical ones.

## **B. Instrumentation Table Preparation Step**

- B.1. Prepare four sterile basins with at least 500 mL of sterile physiological saline in each to rinse the valve prosthesis.
- B.2. Prepare one sterile basin with at least 500 mL of sterile physiological saline for rinsing and flushing the other devices.
- B.3. Prepare a sterile basin with a mixture of 90 mL of sterile physiological saline and 30 mL of contrast solution, to flush the balloon lines of the delivery system and the pre-dilation balloon, as well as to inflate their balloons.
- B.4. Prepare two balloon inflation devices, two three-way stopcocks, a 20 mL syringe, a 60 mL syringe, tweezers and a scalpel.

## **C. Aureus® Transfemoral Expandable Introducer Preparation Step**

- C.1. Remove the Aureus® Transfemoral Expandable Introducer (expandable introducer) from its packaging. The packaging contains an expandable sheath, a dilator and an expander, as shown in Figure 8.

**CAUTION: THE AUREUS® TRANSFEMORAL EXPANDABLE INTRODUCER IS ARRANGED IN A TRAY PACKED IN A STERILE BARRIER ENVELOPE. OPEN THE ENVELOPE AND TRANSFER THE COMPONENTS INSIDE THE TRAY TO THE STERILE FIELD.**

- C.2. Connect a syringe with sterile physiological saline to the expandable sheath's flush port.

3. Inject the solution until it fully fills the inner lumen of the expandable sheath.

4. Remove the syringe from the flush port.

**CAUTION: CLOSING THE FLUSH PORT BEFORE REMOVING THE SYRINGE IS ESSENTIAL TO PREVENT THE LOSS OF STERILE PHYSIOLOGICAL SALINE AND THE ENTRY OF AIR INTO THE SHEATH, WHICH MAY RESULT IN AIR EMBOLISM.**

- C.5. Flush the inner lumen of the dilator and rinse its external surface using sterile physiological saline.

- C.6. Flush the inner lumen of the expander and rinse its external surface using sterile physiological saline.

C.7. Rinse the expandable sheath with sterile physiological saline to hydrate its surface.

**WARNING: THE HYDRATION OF THE EXTERNAL SURFACES AND FILLING OF THE LUMENS OF THE EXPANDABLE SHEATH, DILATADOR AND EXPANDER ARE ESSENTIAL TO INCREASE THEIR HYDROPHILIC PROPERTY, FACILITATING THEIR HANDLING.**

C.8. Insert the dilator completely into the lumen of the expandable sheath and lock the two devices at their hubs, as shown in Figura 17.


Figure 17:  
Locking between the hubs  
of the expandable sheath  
and dilator: a) dilator's hub  
approaching the  
expandable sheath's hub, b)  
hubs connected, needing to  
rotate the dilator hub to  
lock, and  
c) locked hubs.

#### D. INOVARE® SAFESYNC Biological Valvular Prosthesis Preparation Step

D.1. The package containing the INOVARE® SAFESYNC Biological Valvular Prosthesis (valve prosthesis) should not be opened before assuring that the diameter of the device matches the diameter chosen for use in the procedure.

- D.2. Check the temperature indicators on the outer packaging. There are two indicators for recording exposure to extreme temperatures. They should be transparent, indicating that the valve prosthesis was not exposed to temperatures outside the recommended limits during transportation and storage. If any of the indicators are red, then the valve prosthesis has been exposed to temperatures beyond the recommended limits and should not be used; in this case, discard and replace with another valve prosthesis that has not been exposed to temperatures outside the recommended limits.
- D.3. Remove the medical device from the cardboard paper box and styrofoam box, and carefully check for damage or defects, such as cracks in the container or lid, leakage of storage solution, and violation or lack of plastic casing and/or safety seal.
- D.4. To remove the valve prosthesis from the plastic container, peel the plastic casing, remove the safety seal and open the container.
- WARNING: THE OUTER SURFACE OF THE CONTAINER SHOULD NOT BE  
CONSIDERED STERILE. NEVER PLACE IT ON THE INSTRUMENTATION TABLE.**
- D.5. Remove the valve prosthesis from the container using sterile technique, grasping it by the holder with the aid of tweezers.
- CAUTION: DO NOT TOUCH THE BIOLOGICAL TISSUE (CUSPS) OR THE SKIRT FABRIC OF THE VALVE PROSTHESIS WITH THE TWEEZERS, UNDER THE RISK OF CAUSING DAMAGE AND COMPROMISING THE PERFORMANCE OF THE DEVICE.**
- D.6. Cut the three suture lines of the holder assembly with the aid of a scalpel, as illustrated in Figure 18, and completely remove the suture lines to separate the valve prosthesis from its support.


Figure 18: Cutting procedure of the support lines of the valve prosthesis holder assembly.

**CAUTION: DO NOT INSERT THE VALVE PROSTHESIS HOLDER OR THE SUPPORT LINES IN THE SAME RINSING BASIN AS THE VALVE PROSTHESIS; RATHER, DISCARD THEM.**

D.7. After separation of the holder assembly, in the following steps where the valve prosthesis is handled, it should be done exclusively by hand, in order to reduce the risk of trauma and damage to the structures of the valve prosthesis.

**WARNING: THE VALVE PROSTHESIS IS STERILIZED BY GLUTARALDEHYDE AND PRESERVED IN FORMALDEHYDE, NOT REQUIRING ANY TYPE OF ANTIBIOTIC OR FUNGICIDE.**

D.8. The valve prosthesis should be immersed in the first basin containing sterile physiological saline and rinsed for at least two minutes. To rinse, gently and slowly swirl the valve prosthesis back and forth. During rinsing, make sure that the valve prosthesis is fully submerged in the solution.

D.9. After rinsing, transfer the valve prosthesis to the second sterile basin and repeat the rinsing process.

- D.10. Repeat the same procedure in the third and fourth sterile basins.

**WARNING: AFTER RINSING AND REMOVING THE VALVE PROSTHESIS FROM A BASIN, DO NOT RETURN IT TO THAT SAME BASIN AGAIN.**

**CAUTION: DO NOT ALLOW THE VALVE PROSTHESIS TO COME IN CONTACT WITH THE SURFACES OF THE BASINS WHEN RINSING.**

**CAUTION: DO NOT INSERT OTHER DEVICES INTO THE SOLUTION IN WHICH THE RINSING OF THE VALVE PROSTHESIS IS BEING OR HAS BEEN CARRIED OUT.**

- D.11. The valve prosthesis should be left in the last sterile rinsing basin until the time of its use (crimping), in order to keep the tissue of the cusps hydrated.

**CAUTION: DO NOT ALLOW THE VALVE PROSTHESIS TO DRY AT ANY TIME.**

## **E. Pre-Dilation Balloon Catheter Preparation Step**

- E.1. Remove the Pre-Dilation Balloon Catheter (pre-dilation balloon) from its packaging.

**CAUTION: THE PRE-DILATION BALLOON CATHETER IS ARRANGED IN A TRAY PACKED IN A STERILE BARRIER ENVELOPE. OPEN THE ENVELOPE AND TRANSFER THE COMPONENTS INSIDE THE TRAY TO THE STERILE FIELD.**

- E.2. Flush the guidewire line using a syringe filled with sterile physiological saline.

- E.3. Connect a three-way stopcock to the balloon inflation port.

- E.4. Fill the syringe with 40 mL of the mixture of sterile saline and contrast solution and connect it to the three-way stopcock.

- E.5. Connect a balloon inflation device to the three-way stopcock.

- E.6. Close the three-way stopcock to the balloon inflation device port.

- E.7. Using the syringe, fill the balloon catheter line without fully unwrapping the balloon (partial inflation). Retrieve the solution back into the syringe, releasing the plunger in a controlled manner. Repeat the procedure while the presence of air is observed on the line, so as to leave zero pressure in the system.

- E.8. Close the three-way stopcock to the pre-dilation balloon.

E.9. Transfer to the balloon inflation device only the volume necessary to perform the pre-dilation, according to Table 8.

**CAUTION: INSPECT THE INFLATION DEVICE TO ENSURE THAT THE DEVICE IS COMPLETELY FREE OF AIR BUBBLES.**

**CAUTION: USE ONLY THE INDICATED INFLATION MEDIUM. NEVER USE AIR OR OTHER GASEOUS MEDIA TO INFLATE BALLOON DEVICES**

Table 8: Relationship between the diameter of the pre-dilation balloon and the respective nominal volume of inflation.

Pre-dilation balloon catheter	Nominal inflation volume
18 mm	11 mL
20 mm	14 mL

E.10. Close the three-way stopcock to the syringe port;

E.11. Disconnect the syringe from the three-way stopcock.

## F. Crimper Preparation Step

F.1. Remove the crimper from its packaging;

**CAUTION: THE CRIMPER IS PACKED IN A STERILE DOUBLE MICROBIAL BARRIER ENVELOPE. OPEN ONLY THE OUTER ENVELOPE AND TRANSFER THE CRIMPER INSIDE THE INNER ENVELOPE TO THE STERILE FIELD, WHERE IT CAN THEN BE REMOVED FROM THE INNER ENVELOPE.**

F.2. Check the integrity of the compression system by fully turning the handle along its course, in order to observe the correct and total opening of the crimper orifice, that is, its crimping area, as well as the correct and total closure of the crimping area.

- F.3. Rinse the crimper with physiological saline.
- F.4. Position the crimper on the sterile table supported on a stable base (e.g., ends of the table, on top of one of the legs), in order to achieve better stability during crimping of the valve prosthesis.

## **G. B-Delivery System Preparation Step**

- G.1. Remove the B-Delivery System (delivery system) from its packaging.

**CAUTION: THE B-DELIVERY SYSTEM IS PACKED IN A STERILE DOUBLE MICROBIAL BARRIER ENVELOPE. OPEN ONLY THE EXTERNAL ENVELOPE AND TRANSFER THE DELIVERY SYSTEM INSIDE THE INTERNAL ENVELOPE TO THE STERILE FIELD, WHERE THEN THE DELIVERY SYSTEM WILL BE REMOVED FROM THE INTERNAL ENVELOPE.**

- G.2. Check the identification of the diameter of the inflation balloon of the delivery system, located on its handle, ensuring that it is in agreement with the diameter of the INOVARE® SAFESYNC Biological Valvular Prosthesis to be implanted.

- G.3. Remove the protection tube from the delivery system.

- G.4. Remove the lumen protection wire by theatraumatic tip of the delivery system.

- G.5. Flush the delivery system with sterile physiological saline through the flush port, using a syringe.

- G.6. Flush the guidewire line using a syringe, filling with sterile physiological saline through the correct lumen in the hub.

- G.7. Check the integrity of the positioning system:

- G.7.1. Keeping the positioning rod fixed in position with the right hand, hold the handle with the left hand use the left thumb to press the positioning rod unlocking button.

- G.7.2. With the button pressed, start moving the handle towards the hub, releasing the button button just after the movement of the handle starts, and fully pull it back until it touches the hub.

- G.7.3. During this process, monitor and observe the axial displacement of the sheath of the delivery system.

- G.7.4. Initiate the advancement of the handle (moving away from the hub), and, upon reaching the crimping and deployment position, the button will "click", indicating the position and locking the movement of the handle.

- G.7.5. Press the button and once again advance the handle, releasing the button just after the handle movement starts, until another "click" is heard, now indicating the delivery position.
- G.7.6. During these processes, monitor and observe the axial displacement of the sheath of the delivery system.
- G.7.7. Return the delivery system to the crimping and deployment position by pressing the button, starting the retrieval, and releasing the button just after the movement of the handle starts, until the "click" is heard indicating the position.
- G.8. Check the integrity of the deflection system: completely rotate the flex wheel located on the handle, observing the deflectable sheath's deflection. After checking, return the flex wheel to the starting position. Observe the markings near the flex wheel to check the correct orientation to rotate it, as shown in Figure 19.


Figure 19: Flex wheel markings indicating the turning direction to remove deflection (counterclockwise) and apply deflection (clockwise).

- G.9. Check the integrity of the fine adjustment system: completely rotate the fine adjustment wheel located at the base of the handle, observing the axial displacement of the inflation balloon, and then return the fine adjustment wheel to its initial position. Observe the markings near the fine adjustment wheel to check the correct orientation to rotate, as shown in Figure 20.


Figure 20:  
Fine adjustment  
wheel markings  
indicating the turning  
direction to advance  
the balloon  
(counterclockwise)  
and retrieve the  
balloon (clockwise).

**WARNING: IF ANY OF THE ABOVE TESTS ARE CONSIDERED UNSATISFACTORY,  
REPLACE THE DELIVERY SYSTEM AND PERFORM THE STEPS ABOVE FOR THE  
NEW DELIVERY SYSTEM.**

- G.10. Connect a three-way stopcock to the balloon inflation hub.  
G.11. Fill the syringe with 40 mL of the mixture of sterile physiological saline and contrast medium and connect it to the three-way stopcock.  
G.12. Connect a balloon inflation device to the three-way stopcock.  
G.13. Close the three-way stopcock to the balloon inflation device.  
G.14. Using the syringe, fill the inflation balloon line without fully unwrapping the balloon (partial inflation). Retrieve the solution back into the syringe, releasing the plunger in a controlled manner. Repeat the procedure while the presence of air is observed on the line, so as to leave zero pressure in the system.

- G.15. Close the three-way stopcock to the delivery system.
- G.16. Leave in the balloon inflation device only the necessary volume to perform the implantation of the valve prosthesis, according to Table 9.

Table 9: Relationship between the balloon diameter of the B-Delivery System and its nominal inflation volume.

Diameter of the delivery system's inflation balloon	Nominal inflation volume
20 mm	11 mL
22 mm	15 mL
24 mm	21 mL
26 mm	24 mL
28 mm	27 mL
30 mm	30 mL

**CAUTION: INSPECT THE BALLOON INFLATION DEVICE TO ENSURE THAT THE DEVICE (BARREL AND EXTENSION LINE) IS COMPLETELY FREE OF AIR BUBBLES.**

**CAUTION: USE ONLY THE INDICATED INFLATION MEDIUM. NEVER USE GASEOUS OR AIR MEDIA TO INFLATE BALLOON DEVICES.**

- G.17. Close the three-way stopcock to the delivery system.
- G.18. Reinsert the lumen protection wire in the atraumatic tip of the delivery system.

## H. Valve Prosthesis Loading Step

- H.1. Make sure that the delivery system is in the *crimping and deployment position* and that the lumen protection wire is inserted.
- H.2. Position the valve prosthesis between the radiopaque markers of the inflation balloon, **with its aortic face pointing towards the handle of the delivery system**, as shown in Figure 21.


Figure 21: Correct orientation of the INOVARE® SAFESYNC Biological Valvular Prosthesis on the B-Delivery System inflation balloon - with the aortic face pointing towards the delivery system's handle.

- H.3. Insert the delivery system with the valve prosthesis in the crimper's crimping area. Ensure that the valve prosthesis is positioned as close as possible to the center of the crimping area and that the atraumatic tip of the delivery system is completely outside the crimping area.
- H.4. Perform partial crimping of the valve prosthesis, rotating the crimper's handle up to about 75% of the total course, as shown in Figure 22, and then return it to the initial position.


Figure 22: Partial crimping of the INOVARE® SAFESYNC Biological Valvular Prosthesis in the B-Delivery System, turning the Crimper's handle to about 75% of its total course.

- H.5. Keeping the delivery system with the pre-crimped valve prosthesis fixed in position, move the crimper away in a controlled manner, completely removing the system from the crimping area.
- H.6. Adjust the positioning of the valve prosthesis on the inflation balloon. For correct positioning, the valve prosthesis should be about 5 mm from the distal end of theatraumatic tip ("legs" of the nose cone), as shown in Figure 23.


Figure 23: Correct positioning of the INOVARE® SAFESYNC Biological Valvular Prosthesis on the B-Delivery System balloon – at about 5 mm from the distal end of the tip ("legs" of the nose cone).

- H.7. Insert the delivery system with the pre-crimped valve prosthesis in the crimping area. Ensure that the valve prosthesis is positioned as close as possible to the center of the crimping area and that theatraumatic tip of the delivery system is completely outside the crimping area.
- H.8. Perform total crimping of the valve prosthesis, completely rotating the crimper's handle and keeping it in position for fifteen seconds. Return the handle to the starting position.
- H.9. Keeping the delivery system with the crimped valve prosthesis fixed in position, move the crimper away in a controlled manner, completely removing the system from the crimping area.
- H.10. Perform the crimped diameter test, inserting the inflation balloon with the crimped valve prosthesis in the small hole located at the base of the crimper, and verifying whether the crimped valve prosthesis can cross it or not, as shown in Figure 24.


Figure 24: Illustration of the crimped diameter test, indicating the location of the hole in the base of the crimper.

**WARNING: IF THE CRIMPED DIAMETER TEST IS UNSATISFACTORY, CRIMP THE VALVE PROSTHESIS AGAIN AND REPEAT THE TEST (STEPS , , AND ). IF NECESSARY, CRIMP AND TEST ONE LAST TIME. IF THE TEST RESULT IS UNSATISFACTORY AFTER CRIMPING THE VALVE PROSTHESIS THREE TIMES, REPLACE THE DELIVERY SYSTEM AND THE VALVE PROSTHESIS AND PERFORM ALL THE STEPS FOR THE NEW DEVICES.**

- H.11. Perform the negative pressure test of the delivery system with the crimped valve prosthesis, with the aid of the syringe connected to the three-way stopcock.

**WARNING: IF THE NEGATIVE PRESSURE TEST FAILS (NO VACUUM), REPLACE THE VALVE PROSTHESIS AND THE DELIVERY SYSTEM, BECAUSE THE INFLATION BALLOON IS RUPTURED. PERFORM ALL THE ABOVE STEPS FOR THE NEW DELIVERY SYSTEM AND FOR THE NEW VALVE PROSTHESIS.**

- H.12. Close the three-way stopcock to the syringe, and disconnect it from the three-way stopcock.

- H.13. Position the delivery system in the delivery position, fitting the uncovered stent of the valve prosthesis into the tapered section of the deflectable sheath, as shown in Figure 25. If a better fit is required, perform position adjustment with the fine adjustment system.


Figure 25: Delivery configuration of the delivery system, with the uncovered stent of the valve prosthesis loaded in the tapered section of the deflectable sheath: a) balloon not retrieved and uncovered stent not loaded, and b) uncovered stent loaded.

**WARNING: PERFORM THE LOADING SMOOTHLY AND CAREFULLY, ENSURING THAT THERE IS NO MOVEMENT OF THE VALVE PROSTHESIS (PUSHED BY THE DEFLECTABLE SHEATH).**

- H.14. Activate the hydrophilic coating of the deflectable sheath, moisturizing it with a gauze saturated with sterile physiological saline.

**WARNING: IN CASE THE TIME INTERVAL BETWEEN CRIMPING THE VALVE AND BEGINNING THE DELIVERY IS GREATER THAN FIVE MINUTES, MOISTURIZE THE CRIMPED VALVE PROSTHESIS, WRAPPING IT IN A GAUZE MOISTENED WITH STERILE PHYSIOLOGICAL SALINE.**

**CAUTION: IN ORDER TO PREVENT POSSIBLE DAMAGE TO THE CUSPS, THE VALVE PROSTHESIS MUST NOT REMAIN CRIMPED FOR MORE THAN FIFTEEN MINUTES.**

## I. Aureus® Transfemoral Expandable Introducer Vascular Access Step

**WARNING: CHECK THAT THE CONDITION AND DIAMETER OF THE ACCESS VESSEL ARE SUITABLE FOR THE DIAMETER OF THE INTRODUCER BEFORE AND AFTER ITS EXPANSION, AS THE INCORRECT PROPORTION CONFIGURES SERIOUS RISK OF VASCULAR DAMAGE.**

- I.1. If necessary, pre-dilate the femoroiliac vessel.
- I.2. Use the standard Seldinger technique to establish transfemoral access. A 0.035" guide wire with flexible distal tip must cross the puncture needle to achieve transfemoral access.
- I.3. Gently advance the set formed by the expandable sheath and the dilatador as a single piece over the wire until the expandable sheath's hemostatic sleeve is approximately 10 mm inside the blood vessel.
- I.4. Inject the radiopaque contrast medium through the dilator drain port, thus ensuring correct positioning.
- I.5. Slowly and cautiously remove the dilator from the expandable sheath, in order to minimize blood loss.

**CAUTION: IF THERE IS RESISTANCE DURING INSERTION, HALT THE INTRODUCTION OF THE ASSEMBLY FORMED BY THE EXPANDABLE SHEATH AND THE DILATOR UNTIL THE CAUSE OF THE RESISTANCE IS DETERMINED AND CORRECTED. IF RESISTANCE IS STILL FELT, STOP THE PROCEDURE, REMOVE THE ASSEMBLY AND CHANGE THE DEVICE.**

- I.6. Insert the expander in the expandable sheath over the wire until its hub is duly touching the expandable sheath's hub, identically to the process exemplified in Figure 17's steps a) and b).
- I.7. Observe the expander's advancement under fluoroscopic guidance, making sure that the expansion of the sheath is completed;

**CAUTION: IF THERE IS RESISTANCE DURING INSERTION, HALT THE INTRODUCTION OF THE EXPANDER UNTIL THE CAUSE OF THE RESISTANCE IS DETERMINED AND CORRECTED. IF RESISTANCE IS STILL FELT, STOP THE PROCEDURE, REMOVE THE ASSEMBLY AND CHANGE THE DEVICE.**

**CAUTION: DO NOT ATTEMPT TO FOLD A SHEATH THAT HAS ALREADY BEEN EXPANDED.**

**CAUTION: ALWAYS MAINTAIN THE POSITION OF THE EXPANDABLE SHEATH DURING THE PROCEDURE. THE ADVANCE OF AN EXPANDED SHEATH SHOULD BE AVOIDED, AS THIS MOVEMENT MAY NOT BE POSSIBLE AND ANY ATTEMPT WILL POTENTIALLY RESULT IN VASCULAR DAMAGE.**

- I.8. Remove the expander slowly and cautiously, in order to minimize blood loss. It is now possible to introduce devices or catheters through the expanded expandable sheath, always keeping the guidewire central to the devices so that there is no damage to the hemostatic valves.

**CAUTION: MAINTAIN THE POSITION OF THE EXPANDED SHEATH IN THE BLOOD VESSEL DURING ALL OF THE PROCEDURE. DO NOT MOVE IR FORWARD, AS ANY MOVEMENT MAY CAUSE VASCULAR DAMAGE.**

**WARNING: IF THE EXPANDABLE SHEATH IS ONLY PARTIALLY EXPANDED, STOP THE PROCEDURE. FORCING ANY ITEM THROUGH THE PARTIALLY EXPANDED SHEATH MAY RESULT IN DAMAGE SUCH AS RUPERTING OR JAMMING THE DEVICE.**

## J. Pre-dilation Balloon Access And Valvuloplasty Step

- J.1. Insert the guidewire through the tip of the pre-dilation balloon.
- J.2. Advance the pre-dilation balloon through the hemostatic valves of the expandable sheath, accompanying navigation under fluoroscopic guidance, until the balloon completely crosses the expandable sheath.

**CAUTION: THE USE OF EXCESS CONTRAST SOLUTION MAY LEAD TO ACUTE KIDNEY INJURY. MONITOR THE PATIENT'S CREATININE LEVEL BEFORE THE PROCEDURE. ADMINISTRATION SHOULD BE MONITORED THROUGHOUT THE PROCEDURE.**

- J.3. Continue with the navigation of the pre-dilation balloon to the ascending aorta.
- J.4. Cross the stenotic aortic valve in a controlled manner, keeping the guidewire centered in the valve.
- J.5. Position the pre-dilation balloon so that the virtual ring is centered between the distal and proximal radiopaque markers.
- J.6. With the position defined, activate the pacemaker at high frequency and then inflate the pre-dilation balloon in a controlled manner, until the plunger of the inflation device reaches the graduation of 0 mL / cc (injection of all the volume).

**WARNING: IN CASE THE INSUFLATION IS NOT CARRIED OUT IN A CONTROLLED MANNER, AND/OR IF THE NOMINAL INFLATION VOLUME IS NOT FOLLOWED, IT IS POSSIBLE THAT THE BALLOON BURSTS (VOLUME ABOVE THE NOMINAL) AND THAT THERE IS DIFFICULTY IN REMOVING IT, OR THAT THE AORTIC VALVE IS NOT PROPERLY EXPANDED (VOLUME BELOW THE NOMINAL), NEEDING TO REPEAT THE OPERATION.**


- J.7. Deflate the pre-dilation balloon, retrieving all inflation media back to the inflation device.
- J.8. With the balloon properly deflated, turn off the pacemaker.
- J.9. Retrieve the pre-dilation balloon in a controlled manner through the aorta and the expandable sheath, completely removing it.

## **K. Delivery System Vascular Access and Valve Prosthesis Implantation Step**

- K.1. Move the loader carefully until it is close to the atraumatic tip of the delivery system, not surpassing it as to keep it visible for the later introduction of the guidewire. When passing the loader over the crimped valve prosthesis, perform the advance smoothly and carefully, so that there is no movement of the crimped valve prosthesis (pushed by the loader). If necessary, lightly open the loader to allow easy passage over the crimped valve prosthesis.
- K.2. Remove the lumen protection wire by the atraumatic tip of the delivery system.
- K.3. Insert the guidewire through the atraumatic tip of the delivery system.
- K.4. Advance the loader over the delivery system in order to position the valve prosthesis in the rigid part of the loader.
- K.5. Make sure that the "BRAILE" logo on the delivery system handle is positioned upwards (direct viewing).

**CAUTION: MAINTAIN THE CORRECT ORIENTATION OF THE DELIVERY SYSTEM, WITH THE "BRAILE" LOGO ON ITS HANDLE POSITIONED UPWARDS, DURING THE ENTIRE PROCEDURE (VASCULAR ACCESS, NAVIGATION, IMPLANTATION OF THE VALVE PROSTHESIS AND VASCULAR REMOVAL).**

- K.6. Insert the loader into the expandable sheath until its rigid part meets the hub of the expandable sheath;
- K.7. Advance the delivery system through the loader and the hemostatic valves of the expandable sheath, observing under fluoroscopic guidance, until the inflation balloon with the crimped valve prosthesis completely crosses the expandable sheath;
- K.8. Continue with the navigation of the delivery system, and, when approaching the beginning of the descending aorta, simultaneously deflect and advance the delivery system through the aortic arch, with extreme care and control;
- K.9. In the ascending aorta, verify the positioning of the distal radiopaque marker of the inflation balloon, which should overlap with the radiopaque marker of the deflectable sheath. If the balloon advances during navigation, reposition using the delivery system's fine-tuning system, as shown in Figure 26.

**Figure 26:**  
Overlap between the distal radiopaque marker of the inflation balloon and the radiopaque marker of the deflectable sheath.

**WARNING: THE OVERLAP OF THE INFLATION'S BALLOON DISTAL RADIOPAQUE MARKER WITH THE RADIOPAQUE MARKER OF THE DEFLECTABLE SHEATH IS INDISPENSABLE FOR THE SYSTEM TO HAVE THE MAXIMUM COLUMNAR FORCE TO CROSS THE CALCIFIED AORTIC VALVE.**

- K.10. Cross the aortic valve in a controlled manner, keeping the guidewire centered. If necessary, use the deflection system to adjust the coaxiality of the delivery system in the aortic valve.
- K.11. After crossing the valve, return the fine adjustment wheel to its initial position to undo all displacement of the inflation balloon (Note: ignore this step if fine adjustment was not applied);
- K.12. Position the valve prosthesis, aligning its proximal radiopaque marker with the non-coronary valve sinus, as shown in Figure 27.


Figure 27: Alignment between the proximal radiopaque marker of the valve prosthesis and the aortic valve sinus. The right figure presents visual guides for the radiopaque marker and valve sinus, being NCC = non-coronary cusp, RCC = right coronary cusp and LCC = left coronary cusp. The yellow dashed line indicates the alignment plane.

- K.13. Perform the pullback maneuver: keeping the positioning rod fixed, press the unlocking button and retrieve the handle towards the balloon inflation port, releasing the locking button just after starting of the movement, and bring the handle to the deployment position. Upon reaching the position, the button will "click", locking the position rod. This movement allows for the full exposure of the inflation balloon, as shown in Figure 28.


Figure 28: Fluoroscopic images demonstrating the pullback maneuver: a) before and b) after. Note the distance of the deflectable sheath from the inflation balloon (distance between the sheath's marker and the distal marker of the balloon).

- K.14. The positioning of the valve prosthesis in the aortic valve can be improved using the fine adjustment of the delivery system, as needed.
- K.15. If necessary, use the deflection system to adjust the coaxiality of the valve prosthesis in the aortic valve.
- K.16. With the position set, activate the pacemaker at high frequency inflate the balloon in a controlled manner, until the plunger of the inflation device reaches the graduation of 0 mL / 0 cc (injection of the entire volume).

**WARNING: IF THE INFLATION IS NOT CARRIED OUT IN A CONTROLLED MANNER, AND/OR IF THE NOMINAL INFLATION VOLUME IS NOT FOLLOWED, IT IS POSSIBLE THAT:**

- **VOLUME ABOVE THE NOMINAL: THE BALLOON BURSTS, ELEVATING THE RISKS DURING RETRIEVAL, OR THAT THERE IS RUPTURE OF AORTIC RING, PARAVALVULAR LEAK AND/OR NEED FOR PACEMAKER.**
- **VOLUME BELOW THE NOMINAL: THE VALVE PROSTHESIS MAY NOT BE FULLY EXPANDED, LEADING TO MIGRATION OR REQUIRING SUBSEQUENT POST-DILATATION.**

- K.17. Deflate the inflation balloon, retrieving the solution the inflation balloon.
- K.18. With the balloon properly deflated, turn off the pacemaker.

## **L. B-Delivery System Vascular Removal Step**

- L.1. Initiate the retrieval of the delivery system of the aortic valve in a controlled manner, keeping the guidewire central to the valve, in order to avoid possible displacement of the implanted valve prosthesis.
- L.2. With the inflation balloon positioned in the ascending aorta, remove all deflection from the deflectable sheath and return the fine adjustment wheel to its initial position.
- L.3. Retrieve the delivery system through the ascending aorta and aortic arch to the descending aorta.
- L.4. With the inflation balloon positioned in the descending aorta, perform the complete inflation balloon retrieval, keeping the handle fixed in its position, pressing the unlocking button and retreating the positioning rod to the limit of its course. In addition, apply fine adjustment by fully rotating the thread.
- L.5. Perform complete removal of the delivery system by retreating it through the expandable sheath.

## M. Positioning Verification and Measurement of Valve Prosthesis

- M.1. Perform a supraaortic angiography to evaluate device performance and coronary clearance. The ideal position should meet the criteria illustrated in Figure 29.


Figure 29: Ideal positioning of the INOVARE® SAFESYNC Biological Valvular Prosthesis. In the image on the left, it is observed the positioning of 25% of the valve in the ventricular portion, and the remaining 75% in the aortic portion, respecting the maximum required depth. The goal is to provide a mostly supra-annular implant, ensuring better performance and reducing the risk of conduction disorders. In the image on the right there are visual guides of the alignment of the valve prosthesis implanted in relation to the valve sinus.

- M.2. Monitor and record transvalvular pressure gradients and evaluate the performance of the valve prosthesis.  
M.3. After satisfactory distribution, remove all devices when the TCA level is appropriate.  
M.4. Perform an angiography to check for any vascular complications.  
M.5. Close the access location.

## **9. VALIDITY**

The INOVARE® SAFESYNC Biological Valvular Prosthesis - B-Delivery System is valid for 2 years.

Validity is counted from the date of sterilization, and the product is considered sterile in this period as long as package integrity is maintained.

## **10. LIMITED WARRANTY**

BRAILE BIOMÉDICA INDÚSTRIA, COMÉRCIO E REPRESENTAÇÕES LTDA., attests that the product was manufactured with the proper care, with good quality materials, and is guaranteed against manufacturing vices in accordance with the Consumer's Defense Code. It clarifies further that, if the product is tampered by unauthorized or unskilled individuals, the warranty will lose its validity. The use of the product is restricted to medical or qualified and trained technicians under the supervision of the doctor in charge, who are familiar with the procedures inherent to the product, and being responsible for the product's use. The manufacturer shall not be responsible for any damage arising from improper use of the product. BRAILE BIOMÉDICA LTDA establishes further that, the right of complaints with respect to the product will be governed by the dispositions contained in the Consumer's Defense Code (Law 8.078/90).

No agent, employee, representative or distributor of BRAILE BIOMÉDICA LTDA has authority to change or amend what herein is described, assume or link BRAILE BIOMÉDICA LTDA to any other responsibility or warranty in connection to this product.

## 11. DESCRIPTION OF SYMBOLS

	Consult instructions for use		Use-by-date
<b>STERILE</b>	Sterile	<b>REF</b>	Catalogue number
<b>STERILE EO</b>	Sterilized using ethylene oxide		Temperature limit
	Non-pyrogenic		Fragile, handle with care
	Do not re-use		This side up
	Do not use if package is damaged		Keep dry
	Do not resterilize		Keep away from sunlight
<b>SN</b>	Serial number		Stacking limit by number
<b>LOT</b>	Batch code		Recyclable
	MR conditional	<b>MD</b>	Medical device
	Date of manufacture		
	Manufacturer		

## 12. LABEL MODEL – Corrugated carton box containing INOVARE® SAFESYNC Biological Valvular Prosthesis, Pre-Dilation Balloon Catheter, Aureus® Transfemoral Expandable Introducer and Crimper

	Braille Biomédica Indústria Comércio e Representações Ltda. Av. Presidente Juscelino Kubitschek de Oliveira, 1505 - Jardim Tarraf I ZIP Code 15091-450 - São José do Rio Preto - SP - Brazil Phone: +55 17 2136-7000 SAC: 0800 707 20 50 CNPJ: 52.828.936/0001-09 www.braille.com.br		
<b>Technical Name:</b> Biological Cardiac Valve Prosthesis			
<b>Trade Name:</b> INOVARE® SAFESYNC Biological Valvular Prosthesis – B-Delivery System			
<b>Model:</b> xxxxx			
<b>ANVISA registration no.</b> 10159030114			
<b>Technician in Charge:</b> Vladimir D. A. Ramirez - CRF-SP 09010			
<b>REF</b>	<b>SN</b>	<b>LOT</b>	
<b>Packaging contents:</b> 01 (one) INOVARE® SAFESYNC Biological Valvular Prosthesis (20, 22, 24, 26, 28 or 30 mm); 01 (one) Pre-Dilation Balloon Catheter (18 or 20 mm); 01 (one) Aureus® Transfemoral Expandable Introducer; 01 (one) Crimper, 01 (one) Instruction Manual.			
Single-use product. Do not reprocess.			
Liquid chemical sterilization (valve prosthesis) or sterilized by ethylene oxide (other components).			
	DD/MM/YY		DD/MM/YY
	STERILE		STERILE
			
	Warnings and precautions, technical characteristics and conditions for storage, transportation and handling the product: consult instructions for use.		
The product, in its original package, should be stored in a clean, dry, weather-free environment and at temperatures above 10 °C and below 30 °C. Storage also requires proper care to prevent falls and possible damage to the product. For transport, it should never be subjected to freezing or exposed to extreme temperatures. During transportation, avoid exposure to water, direct sunlight, extreme temperatures and high humidity.			
			
			
For the internationally harmonized meanings for symbols - consult instructions for use.			

## 13. LABEL MODEL – Corrugated carton box containing B-Delivery System

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**Technical Name:** Biological Cardiac Valve Prosthesis

**Trade Name:** INOVARE® SAFESYNC Biological Valvular Prosthesis – B-Delivery System

**Component:** INOVARE® SAFESYNC Biological Valvular Prosthesis (20, 22, 24, 26, 28 or 30 mm)

**Model:** xxxxx

**ANVISA registration no.** 10159030114

**Technician in Charge:** Vladimir D. A. Ramirez - CRF-SP 09010

**REF**

**SN**

**LOT**

**Packaging contents:** 01 (one) INOVARE® SAFESYNC Biological Valvular Prosthesis (20, 22, 24, 26, 28 or 30 mm); 01 (one) Instruction Manual, 01 (one) Implant Registration form, 06 (six) traceability labels.

**Composition:** transcatheter biological cardiac valve prosthesis – bovine pericardium tissue preserved in formaldehyde solution, cobalt-chromium alloy, polyester, tantalum.

Single-use product. Do not reprocess.

Liquid chemical sterilization.

DD/MM/YY

DD/MM/YY

STERILE

Warnings and precautions, technical characteristics and conditions for storage, transportation and handling the product: consult instructions for use.

The product, in its original package, should be stored in a clean, dry, weather-free environment and at temperatures above 10 °C and below 30 °C. Storage also requires proper care to prevent falls and possible damage to the product. For transport, it should never be subjected to freezing or exposed to extreme temperatures. During transportation, avoid exposure to water, direct sunlight, extreme temperatures and high humidity.

*For the internationally harmonized meanings for symbols: consult instructions for use.*

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**Technical Name:** Biological Cardiac Valve Prosthesis

**Trade Name:** INOVARE® SAFESYNC Biological Valvular Prosthesis – B-Delivery System

**Component:** B-Delivery System (20, 22, 24, 26, 28 or 30 mm)

**Model:** xxxxx

**ANVISA registration no.** 10159030114

**Technician in Charge:** Vladimir D. A. Ramirez - CRF-SP 09010

**REF**

**SN**

**LOT**

**Packaging contents:** 01 (one) B-Delivery System (20, 22, 24, 26, 28 or 30 mm, corresponding to the valve prosthesis model it accompanies).

Composition: delivery system for balloon expandable transcatheter valve prosthesis – polyether-block-amide, polyamide, stainless steel.

Single-use product. Do not reprocess.

DD/MM/YY

DD/MM/YY

**STERILE**

**STERILE EO**

Warnings and precautions, technical characteristics and conditions for storage, transportation and handling the product: consult instructions for use.

The product, in its original package, should be stored in a clean, dry, weather-free environment and at temperatures above 10 °C and below 30 °C. Storage also requires proper care to prevent falls and possible damage to the product. For transport, it should never be subjected to freezing or exposed to extreme temperatures. During transportation, avoid exposure to water, direct sunlight, extreme temperatures and high humidity.

*For the internationally harmonized meanings for symbols: consult instructions for use.*

## 16. LABEL MODEL – Pre-Dilation Balloon Catheter

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**Technical Name:** Biological Cardiac Valve Prosthesis

**Trade Name:** INOVARE® SAFESYNC Biological Valvular Prosthesis – B-Delivery System

**Component:** Pre-Dilation Balloon Catheter (18 or 20 mm)

**Model:** XXXXX

**ANVISA registration no.** 10159030114

**Technician in Charge:** Vladimir D. A. Ramirez - CRF-SP 09010

**Packaging contents:** 01 (one) Pre-Dilation Balloon Catheter (18 or 20 mm).

**Composition:** pre-dilation balloon catheter – polyether-block-amide, polyamide.

Single-use product. Do not reprocess.

DD/MM/YY

DD/MM/YY

Warnings and precautions, technical characteristics and conditions for storage, transportation and handling the product: consult instructions for use.

The product, in its original package, should be stored in a clean, dry, weather -free environment and at temperatures above 10 °C and below 30 °C. Storage also requires proper care to prevent falls and possible damage to the product. For transport, it should never be subjected to freezing or exposed to extreme temperatures. During transportation, avoid exposure to water, direct sunlight, extreme temperatures and high humidity.

For the internationally harmonized meanings for symbols: consult instructions for use.

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**Technical Name:** Biological Cardiac Valve Prosthesis

**Trade Name:** INOVARE® SAFESYNC Biological Valvular Prosthesis – B-Delivery System

**Component:** Aureus® Transfemoral Expandable Introducer

**Model:** xxxxx

**ANVISA registration no.** 10159030114

**Technician in Charge:** Vladimir D. A. Ramirez - CRF-SP 09010

**Packaging contents:** 01 (one) Aureus® Transfemoral Expandable Introducer.

**Composition:** expandable vascular introducer – polyether-block-amide, polytetrafluoroethylene, polyethylene.

**Single-use product.** Do not reprocess.

DD/MM/YY

DD/MM/YY

Warnings and precautions, technical characteristics and conditions for storage, transportation and handling the product: consult *instructions for use*.

The product, in its original package, should be stored in a clean, dry, weather -free environment and at temperatures above 10 °C and below 30 °C. Storage also requires proper care to prevent falls and possible damage to the product. For transport, it should never be subjected to freezing or exposed to extreme temperatures. During transportation, avoid exposure to water, direct sunlight, extreme temperatures and high humidity.

For the internationally harmonized meanings for symbols: consult *instructions for use*.

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**Technical Name:** Biological Cardiac Valve Prosthesis

**Trade Name:** INOVARE® SAFESYNC Biological Valvular Prosthesis – B-Delivery System

**Component:** Crimper

**Model:** xxxxx

**ANVISA registration no.** 10159030114

**Technician in Charge:** Vladimir D. A. Ramirez - CRF-SP 09010

**REF**

**SN**

**LOT**

**Packaging contents:** 01 (one) Crimper.

Composition: transcatheter valve prosthesis crimping device - polyoxymethylene, aluminum alloy.

Single-use product. Do not reprocess.

DD/MM/YY

DD/MM/YY

**STERILE**

**STERILE EO**

**MD**

Warnings and precautions, technical characteristics and conditions for storage, transportation and handling the product: consult *instructions for use*.

The product, in its original package, should be stored in a clean, dry, weather -free environment and at temperatures above 10 °C and below 30 °C. Storage also requires proper care to prevent falls and possible damage to the product. For transport, it should never be subjected to freezing or exposed to extreme temperatures. During transportation, avoid exposure to water, direct sunlight, extreme temperatures and high humidity.

*For the internationally harmonized meanings for symbols: consult instructions for use.*

## 19. LABEL MODEL – INOVARE® SAFESYNC Biological Valvular Prosthesis (Replacement Component)

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**Technical Name:** Biological Cardiac Valve Prosthesis

**Trade Name:** INOVARE® SAFESYNC Biological Valvular Prosthesis – B-Delivery System

**Component:** INOVARE® SAFESYNC Biological Valvular Prosthesis (20, 22, 24, 26, 28 or 30 mm)  
Replacement component for exclusive use in the INOVARE® SAFESYNC Biological Valvular Prosthesis – B-Delivery System.

**Model:** xxxx

**ANVISA registration no.** 10159030114

**Technician in Charge:** Vladimir D. A. Ramirez - CRF-SP 09010

**REF**

**SN**

**LOT**

**Packaging contents:** 01 (one) INOVARE® SAFESYNC Biological Valvular Prosthesis (20, 22, 24, 26, 28 or 30 mm); 01 (one) Instructions Manual, 01 (one) Implant Registration Form, 06 (six) traceability labels.

**Composition:** transcatheter biological cardiac valve prostheses – bovine pericardium tissue preserved in formaldehyde solution, cobalt-chromium alloy, polyester, tantalum.

**Single-use product.** Do not reprocess.

Liquid chemical sterilization.

DD/MM/YY

DD/MM/YY

**STERILE**

Warnings and precautions, technical characteristics and conditions for storage, transportation and handling the product: consult instructions for use.

The product, in its original package, should be stored in a clean, dry, weather-free environment and at temperatures above 10 °C and below 30 °C. Storage also requires proper care to prevent falls and possible damage to the product. For transport, it should never be subjected to freezing or exposed to extreme temperatures. During transportation, avoid exposure to water, direct sunlight, extreme temperatures and high humidity.

For the internationally harmonized meanings for symbols: consult instructions for use.

ENGLISH

## 20. LABEL MODEL – B-Delivery System (Replacement Component)

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**Technical Name:** Biological Cardiac Valve Prosthesis

**Trade Name:** INOVARE® SAFESYNC Biological Valvular Prosthesis – B-Delivery System

**Component:** B-Delivery System (20, 22, 24, 26, 28 or 30 mm)

Replacement component for exclusive use in the INOVARE® SAFESYNC Biological Valvular Prosthesis – B-Delivery System.

**Model:** xxxx

**ANVISA registration no.** 10159030114

**Technician in Charge:** Vladimir D. A. Ramirez - CRF-SP 09010

**Packaging contents:** 01 (one) B-Delivery System (20, 22, 24, 26, 28 or 30 mm, corresponding to the valve prosthesis model it accompanies), 01 (one) Instructions Manual.

**Composition:** delivery system for balloon expandable transcatheter valve prosthesis – polyether-block-amide, polyamide, stainless steel.

Single-use product. Do not reprocess.

DD/MM/YY

DD/MM/YY

Warnings and precautions, technical characteristics and conditions for storage, transportation and handling the product: consult instructions for use.

The product, in its original package, should be stored in a clean, dry, weather-free environment and at temperatures above 10 °C and below 30 °C. Storage also requires proper care to prevent falls and possible damage to the product. For transport, it should never be subjected to freezing or exposed to extreme temperatures. During transportation, avoid exposure to water, direct sunlight, extreme temperatures and high humidity.

*For the internationally harmonized meanings for symbols: consult instructions for use.*

## 21. LABEL MODEL – Pre-Dilation Balloon Catheter (Replacement Component)

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**Technical Name:** Biological Cardiac Valve Prosthesis

**Trade Name:** INOVARE® SAFESYNC Biological Valvular Prosthesis – B-Delivery System

**Component:** Pre-Dilation Balloon Catheter (18 or 20 mm)

Replacement component for exclusive use in the INOVARE® SAFESYNC Biological Valvular Prosthesis – B-Delivery System.

**Model:** xxxxx

**ANVISA registration no.** 10159030114

**Technician in Charge:** Vladimir D. A. Ramirez - CRF-SP 09010

**REF**

**SN**

**LOT**

**Packaging contents:** 01 (one) Pre-Dilation Balloon Catheter (18 or 20 mm), 01 (one) Instructions Manual.

Composition: pre-dilation balloon catheter – polyether-block-amide, polyamide.

Single-use product. Do not reprocess.

DD/MM/YY

DD/MM/YY

**STERILE**

Warnings and precautions, technical characteristics and conditions for storage, transportation and handling the product: consult instructions for use.

The product, in its original package, should be stored in a clean, dry, weather-free environment and at temperatures above 10 °C and below 30 °C. Storage also requires proper care to prevent falls and possible damage to the product. For transport, it should never be subjected to freezing or exposed to extreme temperatures. During transportation, avoid exposure to water, direct sunlight, extreme temperatures and high humidity.

For the internationally harmonized meanings for symbols: consult instructions for use.

## 22. LABEL MODEL – Aureus® Transfemoral Expandable Introducer (Replacement Component)

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Phone: +55 17 2136-7000 SAC: 0800 707 20 50

CNPJ: 52.828.936/0001-09 www.braille.com.br

**Technical Name:** Biological Cardiac Valve Prosthesis

**Trade Name:** INOVARE® SAFESYNC Biological Valvular Prosthesis – B-Delivery System

**Component:** Aureus® Transfemoral Expandable Introducer

Replacement component for exclusive use in the INOVARE® SAFESYNC Biological Valvular Prosthesis – B-Delivery System.

**Model:** xxxxx

**ANVISA registration no.** 10159030114

**Technician in Charge:** Vladimir D. A. Ramirez - CRF-SP 09010

**REF**

**SN**

**LOT**

**Packaging contents:** 01 (one) Aureus® Transfemoral Expandable Introducer, 01 (one) Instructions Manual.

Composition: expandable vascular introducer – polyether-block-amide, polytetrafluoroethylene, polyethylene.

Single-use product. Do not reprocess.

DD/MM/YY

DD/MM/YY

**STERILE**

**STERILE EO**

Warnings and precautions, technical characteristics and conditions for storage, transportation and handling the product: consult *instructions for use*.

The product, in its original package, should be stored in a clean, dry, weather -free environment and at temperatures above 10 °C and below 30 °C. Storage also requires proper care to prevent falls and possible damage to the product. For transport, it should never be subjected to freezing or exposed to extreme temperatures. During transportation, avoid exposure to water, direct sunlight, extreme temperatures and high humidity.

*For the internationally harmonized meanings for symbols: consult instructions for use.*

## 23. LABEL MODEL – Crimper (Replacement Component)

 <p><b>Braile Biomédica Indústria Comércio e Representações Ltda.</b>          Av. Presidente Juscelino Kubitschek de Oliveira, 1505 - Jardim Tarraf I          ZIP Code 15091-450 - São José do Rio Preto - SP - Brazil          Phone: +55 17 2136-7000 SAC: 0800 707 20 50          CNPJ: 52.828.936/0001-09 www.braile.com.br</p>	
<p><b>Technical Name:</b> Biological Cardiac Valve Prosthesis</p> <p><b>Trade Name:</b> INOVARE® SAFESYNC Biological Valvular Prosthesis – B-Delivery System</p> <p><b>Component:</b> Cimper</p> <p>Replacement component for exclusive use in the INOVARE® SAFESYNC Biological Valvular Prosthesis – B-Delivery System.</p> <p><b>Model:</b> XXXXX</p> <p><b>ANVISA registration no.</b> 10159030114</p> <p><b>Technician in Charge:</b> Vladimir D. A. Ramirez - CRF-SP 09010</p>	
<b>REF</b> <b>SN</b> <b>LOT</b>	
<p><b>Packaging contents:</b> 01 (one) Cimper, 01 (one) Instructions Manual.</p> <p><b>Composition:</b> transcatheter valve prosthesis crimping device - polyoxymethylene, aluminum alloy.</p> <p><b>Single-use product. Do not reprocess.</b></p>	
 DD/MM/YY   DD/MM/YY	
     	
 <p>Warnings and precautions, technical characteristics and conditions for storage, transportation and handling the product: consult instructions for use.</p>	
<p>The product, in its original package, should be stored in a clean, dry, weather-free environment and at temperatures above 10 °C and below 30 °C. Storage also requires proper care to prevent falls and possible damage to the product. For transport, it should never be subjected to freezing or exposed to extreme temperatures. During transportation, avoid exposure to water, direct sunlight, extreme temperatures and high humidity.</p>	
       	
<i>For the internationally harmonized meanings for symbols: consult instructions for use.</i>	

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